


## SPECIFICATION SHEET

VERSION 1.1 | AUGUST 2020



# IONIC<sup>®</sup> FLUIDIC CHIP



Ionic<sup>®</sup> Fluidic Chip

### Dimensions


Height	14.53 mm / 0.57 in
Width	85.48 mm / 3.37 in
Length	127.76 mm / 5.03 in
Dispense depth	3 mm below surface

### Working Volumes

Sample well volume	200 $\mu$ L (typical)
Buffer well volume	Varies by protocol, sample type
Elution well volume	40 $\mu$ L (typical)

### Storage Handling

Temperature	15 – 32°C
Packaging	In sealed, static dissipative, bag
Handling	With disposable gloves, avoiding contact with surface



### PURIGEN BIOSYSTEMS, INC.

5700 Stoneridge Drive, Suite 100  
Pleasanton, CA 94588  
USA

TEL: +1 925 264-1364

### SALES

EMAIL: [sales@purigenbio.com](mailto:sales@purigenbio.com)

### SUPPORT

EMAIL: [support@purigenbio.com](mailto:support@purigenbio.com)

TEL: +1 877 787-4436

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2020 Purigen Biosystems Inc. All rights reserved.

Ionic<sup>™</sup> is a trademark of Purigen Biosystems, Inc., in the U.S.

**PURIGEN**<sup>™</sup>  
BIOSYSTEMS

PUR-DOC-15